森佰
SENSBY

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
The instrument SPR-T200 is based on Surface Plasmon Resonance (SPR) technology to monitor the interaction of biomolecules in real time. The molecular interactions information is obtained by measuring the angular reflection spectroscopy at a fixed wavelength; or measuring the kinetic of molecular interactions through reflectivity intensity changes at a fixed incident angle. SPR-T200 is able to perform the determination of kinetic constant, detection of sample concentration; analysis of the interaction modes.

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- ✓ 蛋白质组学研究 (Proteomics)
- ✓ 免疫测定法 (Immunoassay)
- ✓ 新药研发 (Drug Discovery)
- ✓ 疫苗开发 (Vaccine Development)
- ✓ 癌症研究 (Cancer Research)
- ✓ 瞬时结合 (Transient Binding)
- ✓ 信号传递 (Cell Signaling)
- ✓ 配体垂钓 (Ligand Fishing)
- ✓ 多分子复合物的结构和组装 (Multi-molecular Complexes)
- ✓ 结合特异性 (Binding Specificity)
- ✓ 分子识别 (Molecular Recognition)
- ✓ 结构与功能的关系 (Structure-function Relationship)
- ✓ 免疫调节 (Immune Regulation)
- ✓ 酶反应 (Enzyme Reaction)

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As shown in the figure, the laser at the wavelength of 632.8 nm passes through a polarizer to become a transverse magnetically polarized beam and is coupled to glass prism to excite surface plasmon resonance of sensor chip. The reflected light intensity R is detected in real time by the changes of incident light angle θ . The detection of molecular interactions can be achieved by the changes of ΔR or $\Delta \theta$; the molecular interaction kinetic constants can also be obtained by measuring the reflectivity changes at a fixed angle.

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